

Technical data			Document references			
			3EST000202-9005 Sammanställningsritning			
Item ref Item ref desig	Quantity Unit	Identity number	Title Context name	Material, type Data, dimensions	Mass Unit	Remarks
		3EST 64-400	Uppsamlingstank BFM73	400 liter		
1	2.000 ea	3EST 60-12	Nivågivare	Evac, Det.nr. 4000190		Reservdel
2	1.000 ea	3EST 60-13	Värmematta uppsaml.tank	BFM Evac. Det.nr. 4000186		Reservdel
3	1.000 ea	3EST 60-14	Termostat	Evac, Det.nr. 1000176		Reservdel
4	1.000 ea	3EST 59-999	Vakuumenhet BFM	Evac, Det.nr. 5400355		Reservdel
5	0.000 ea	3EST 60-1	Ejektor	Evac, Det.nr. 5406125		Reservdel
6	0.000 ea	3EST 60-3	Filterinsats	Evac, Det.nr. 5475206		Reservdel
7	0.000 ea	3EST 60-4	Vakuumbrytare	Evac, Det.nr. 5716183		
8	1.000 ea	3EST 60-5	Värmematta inkl.termostat	Evac, Det.nr. 4000183		Reservdel
9	1.000 ea	3EST 60-6	Vakuumkanponenttavla BF	Evac, Det.nr. 5472194		Reservdel
10	0.000 ea	3EST000201-1802	Filterregulator	SMC EVHS4000-F04	kg	
11	0.000 ea	3EST000201-1805	Manometer med tryckvakt	GP46-10-02	kg	

Sequence Index: 1

Replaced by:		Replaces:		Stycklista	
				Uppsamlingstank BFM73	
				Released	
				3EST 64-400	
				Revision: <u> A </u> Lang: <u> SV </u> Page: <u> 1/2 </u>	

Prepared
2001-06-13
A Eriksson

Checked

Approved
2001-06-27
T Kunze


ADtranz
 SETRA/SETI

Item ref Item ref desig	Quantity Unit	Identity number	Title Context name	Material, type Data, dimensions	Mass Unit	Remarks
12	0.000 ea	3EST000201-1801	Avstängningsventil	SMC EVHS4000-F04	kg	
13	0.000 ea	3EST 60-9	Magnetventil	Evac, Det.nr. 5472758		Reservdel
14	0.000 ea	3EST 61-174	Pneumatikdetaljer BFM	Evac. Det.nr. 4000905		Reservdel
15	1.000 ea	3EST 72-163	Gummipackning	Evac. Det.nr. 5472197		Reservdel

Administrative data

Data, Dimension 400 liter
Classification Decimal reg 61130

Extended data

Supplier: Evac AB

Suppliers identity;5472191

Sequence Index: 1

Replaced by:		Replaces:		Stycklista		
				 ADtranz	Uppsamlingstank BFM73	Released
		Prepared	2001-06-13		A Eriksson	3EST 64-400
		Checked				
		Approved	2001-06-27	T Kunze	SETRA/SETI	Revision: A